


NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Sheet 1/2 2995847

VERSO CROSS				Plan No. E00055449/E_1.1
Scale: 7:20	Drafted	Verified	Confirmed	
Date	16/03/2023	13/06/2024	13/06/2024	
Initials	SD	SD	PJ	

BOM			
No.	Plan	Article	Qty
1	P00023336	2639354 (39)	1
	P00023305	2639554 (42)	1
	P00023279	2638154 (44)	1
	P00023307	2639754 (46)	1
	P00023309	2639954 (48)	1
	P00023311	2640154 (50)	1
2	P00023337	2639454 (39)	1
	P00023306	2639654 (42)	1
	P00023280	2638254 (44)	1
	P00023308	2639854 (46)	1
	P00023310	2640054 (48)	1
	P00023312	2640254 (50)	1
3	P00052559	2085671 (std.)	2
		2085681 (LB)	2
		20856C5 (LG)	2
		20856A3 (OR)	2
4	X002133	1992424	4
5	X003692	1998763	4
6	X004452	1904058	4

BOM			
No.	Plan	Article	Qty
7	X004606	1905091 (39/44/48)	2
	X004607	1905092 (42/46/50)	2
8	X004620	1905090	2
9	M4 DIN985	1907202	4
10	M8 DIN126	1996214	2
11	M4x16	2990185	4
12	M8 DIN985	1996213	1
13	M6x45 DIN912	2994058	4
14	M6 DIN985	2994063	6
15	3.9x9.5 DIN7981	2992349	4
16	M6x50 DIN912	2995258	2
17	M8x110 DIN912	2995259	1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Sheet 1/2

2995847

VERSO CROSS

Plan No.

E00055449/E_1.1

Scale: 7:20	Drafted	Verified	Confirmed
Date	16/03/2023	13/06/2024	13/06/2024
Initials	SD	SD	PJ

