

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Sheet 1/2

2995860

VERSO FRONT SIDE FRAME B02-B03

Plan No.

E00057375/D1.1

Scale:	1:4	Drafted	Verified	Confirmed
Date	06/04/2023	27/07/2023	27/07/2023	
Initials	SD	SD	PJ	



Detail A


BOM			
No.	Plan	Article	Qty
1	P00023273	2638371 (WH)	1
		2638376 (BL)	1
2	P00023284	2638471 (WH-Left)	1
		2638476 (BL-Left)	1
	P00023286	2638571 (WH-Right)	1
		2638576 (BL-Right)	1
3	P00027229	2643800	1
4	Y016062	1736054	1
5	Y016063	1736154	1
6	Y017422	1747971 (WH)	1
		1747976 (BL)	1
7	X003204	1585100	1
8	X004308	1904129	2
9	X004309	1904133	1
10	X004415	17193	1
11	X004420	1904134	1
12	X004470	1713654	1
13	X004471	1713399	1
14	X004487	1904581	1
15	X004597	1905121	1
16	X004689	1734600	1
17	X004758	1906004	1
18	X004876	1747754	1
19	X004818	1906184	1
20	D3x16 DIN1481	1003791	1

BOM			
No.	Plan	Article	Qty.
21	M6x16	1906622	1
22	---	1908236	1
23	M6x6 DIN914	1909924	3
24	---	1991857	2
25	M6x45 DIN912	2994058	4
26	M6 DIN986	2994059	13
27	M6x40 DIN912	2994317	3
28	M6 DIN125	2994414	7
29	M6x35 DIN7984	2994654	1
30	M6x50 DIN912	2995258	1
31	M6 DIN7349	2995261	1
32	M6x40 DIN7984	2995263	2
33	M8 DIN986	2995398	1
34	M6x25 DIN912	2996887	2
35	M8x40 DIN912	2994915	1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Sheet 2/2

2995860

VERS0 FRONT SIDE FRAME B02-B03				Plan No.
				E00057375/D1.1
Scale: 1:4	Drafted	Verified	Confirmed	
Date	06/04/2023	27/07/2023	27/07/2023	
Initials	SD	SD	PJ	