


2991042

V300 DRIVE SYSTEM HEM2				Plan No. E000652 / F
Sheet 1/2				
Scale: 1:6	Drafted	Verified	Confirmed	
Date	29/09/2010	21/03/2019	21/03/2019	
Initials	DD	SD	PJ	


E000652 / F

Material			
No.	Plan	Article	Quantity
1	X002976	1995343	2
2A	Y021225	2012800	1
2B	Y022137	2025100	1
3	X004810	1742000	4
4	Y017195	1742654 (L)	1
5	Y017206	1742754 (R)	1
6	M6x35 ISO7380	1006057	4
7	M6 DIN986	1008492	4
8	---	1901208	1
9	---	1901745	1
10	61904-2RS1	1994568	4
11	Ø20 DIN471	1995259	4
12A	Y021227	2012900	1
12B	Y021231	2013000	1
12C	Y022136	2025200	1
13A	X005892	G103100	1
13B	X005891	G103000	1
14	X005893	1915848	1
16	M6 DIN7349	1003321	1
17	M6x35	2990215	1

2A Seat width 38 to 52
 2B Seat widths 53 to 60
 12A Seat widths 38 to 44
 12B Seat widths 46 to 52
 12C Seat widths 53 to 60
 13A Seat widths 38 to 44
 13B Seat widths 46 to 60

2991042

E000652 / F

V300 DRIVE SYSTEM HEM2				Plan No.
Sheet 2/2				E000652 / F
Scale: 1:6	Drafted	Verified	Confirmed	
Date	29/09/2010	21/03/2019	21/03/2019	
Initials	DD	SD	PJ	