



E001055 / G

B52 UNIVERSAL SERUM CARRIER				Plan No.
Sheet 1/2				E001055 / G
Scale: 1:3	Drafted	Verified	Confirmed	
Date	16/09/2016	05/05/2017	05/05/2017	
Initials	TT	TT	PJ	

Wheelchair overview			
Model	Clamb type	Hand grip extension	Hand grip
D100	8A - 9A	12A	14C
D200	8A - 9A	12A	14C
D400	8B - 9B	12B	14B
D100 / 30°	8A - 9A	On B90	---
D200 / 30°	8A - 9A	On B90	---
D400 / 30°	8A - 9A	On B90	---
V200	8A - 9A	12C	14A
V300	8A - 9A	12C	14A
V500	8A - 9A	12C	14D
V100D	8A - 9A	12A	14A
V300 / 30°	8A - 9A	On B90	---
V200D	8B - 9B	12B	14B
V300D	8B - 9B	12B	14B
V300DL	8B - 9B	12B	14B
V100	8B - 9B	12B	14B
Eclips+	8A - 9A	12C	14A
Eclips X4	8A - 9A	12C	14A
Eclips X4 Act	8A - 9A	12C	14A
XXL	8A - 9A	12C	14A
Eclips+ 30°	8A - 9A	On B90	---
Eclips X4 30°	8A - 9A	On B90	---
Eclips X4 90°	8A - 9A	On B90	---
Eclips+ D	8A - 9A	12A	14A
708D	8A - 9A	12A	14A
Jazz+	8B - 9B	12B	14A
Jazz S50	8B - 9B	12B	14E
Inovys	8A - 9A	On lateral tube	---
EZ-1	8A - 9A	On lateral tube	---
Navix	---	Not available	---
Forest 3 L44	---	Not available	---
Timix L44	---	Not available	---
Forest 2	8A - 9A	On lateral tube	---
Tracer	8A - 9A	On lateral tube	---
Forest 3 L34 TUB	8A - 9A	On back frame	---
Forest kids	8A - 9A	On back frame	---
Timix L34 TUB	8A - 9A	On back frame	---
Springer	8A - 9A	On back frame	---
Forest 3 L34	8A - 9A	12A	14A
Timix L34	8A - 9A	12A	14A
Forest kids D	8A - 9A	12A	14A
Express VB	8A - 9A	12A	14A
Sqoud	8A - 9A	12A	14A
Rapido	8A - 9A	12A	14A

Material			
No.	Plan	Article	Quantity
1	M6 (Φ20xΦ6.5x1.6)	1997328	1
2	M6x20 DIN912/8.8	1001493	3
3	M5x10	1007063	1
4	M6x10	1900641	1
5	M8x30 DIN912	1992581	1
6	Φ25 M10x1	1999935	1
7	X005631	2036400	1
8A	XP01417	2706799 (Φ22, Φ22,2)	1
8B	X005632	2037299 (Φ25,4)	1
9A	XP01418	2706899 (Φ22, Φ22,2)	1
9B	X005633	2037399 (Φ25,4)	1
10	Y022178	2036554	1
11	Y022180	2036601	1
12A	Y018487	17711	1
12B	Y017202	17421	1
12C	Y018538	17710	1
13	Φ5x18	1001589	1
14A	---	1992600 (Φ22)	1
14B	---	1906244 (Φ25.4)	1
14C	---	3990041 (Φ22.2)	1
14D	Y022063	F007241 (Φ22)	1
14E	---	1908236 (Φ25.4)	1

E001055 / G

<b>B52 UNIVERSAL SERUM CARRIER</b> Sheet 2/2				Plan No.
				E001055 / G
Scale: 1:3	Drafted	Verified	Confirmed	
Date	16/09/2016	05/05/2017	05/05/2017	
Initials	TT	TT	PJ	