


Detail A
Scale: 2:7


EP00457 / D

TURIOS ELECTRICAL BACKREST AND SEAT				Plan No.
				EP00457 / D
Scale: 1:7	Drafted	Verified	Confirmed	
Date	26/11/2019	30/04/2020	30/04/2020	
Initials	DB	MK	TK	

Material			
No.	Plan	Article	Quantity
1	YP07411 (S)	2776954	1
	YP07412 (M)	2777054	1
	YP07413 (L)	2777154	1
2	YP07414 (S)	2777254	1
	YP07415 (M)	2777354	1
	YP07416 (L)	2777454	1
3	YP07417 (S)	2777554	1
	YP07418 (M)	2777654	1
	YP07419 (L)	2777754	1
4	YP07421 (M)	2777954	1
5	YP07682 (S)	2600854	1
	YP07683 (M)	2600954	1
	YP07684 (L)	2601054	1
6	YP07706 (S)	2779754	1
	YP07707 (M)	2779854	1
	YP07708 (L)	2779954	1
7	YP07709 (S)	2780054	1
	YP07710 (M)	2780154	1
	YP07711 (L)	2780254	1
8	LA23	1906790	1
9	LINAK_MOTOR	2990303	1
10	LINX 60A	1912294	1
11	X004942	1907377	2

Material			
No.	Plan	Article	Quantity
12	X004943	1907376	2
13	X004960	1907430	2
14	---	1910745	2
15	1007198	1007198	2
16	1903983	1903983	2
17	Led Head Light	1907340	2
18	Indicator Lamp	1991679	2
19	M10x55 ISO 7380	1998518	1
20	M10 DIN125	1001558	9
21	M10 DIN985	1001683	4
22	M10x45	1908848	3
23	M5x40 DIN933	1001427	2
24	M8 DIN125 BLACK	2994415	16
25	M8 DIN985	1001682	10
26	M8x40	1008739	8
27	M6x45	1991017	4
28	M6 DIN985	2994063	15
29	M6X35_BLACK	2994654	8
30	M6 DIN125	1001564	10
31	M6x16 DIN912	1001492	2
32	Y012123	1675800	2
33	YP08148	2608354	1
34	YP08149	2608454	1
35	M8x16	1991858	2
36	XP02856	GP21400	2
37	6x40	1006111	3

EP00457 / D

TURIOS ELECTRICAL BACKREST AND SEAT				Plan No.
				EP00457 / D
Scale: 1:7	Drafted	Verified	Confirmed	
Date	26/11/2019	30/04/2020	30/04/2020	
Initials	DB	MK	TK	