



Sheet 1/2

LEG REST BZ7 TYPE 708D			Plan No. / Plannr.
BEENSTEUN BZ7 TYPE 708D			E000387/L
Scale / Schaal: 3:10	Drafted / Getekend	Verified / Nagezien	Confirmed / Bevestigd
Date / Datum	13/07/2007	11/09/2017	11/09/2017
Init. / Init.	SD	SD	PJ



Material / Materiaal							
No. / Nr.	Plan	Article / Artikel	Number / Aantal	No.	Plan	Article / Artikel	Number
1	Y000593	F000193	1	21	Large L	1627004	1
2	Y002014	1384200	1		Large R	1627104	1
3	Y007151	1721900	1		Small L	1627204	1
4	Y002209	1389600	1		Small R	1627304	1
5	Y002218	16505..	1		Medium L	1627404	1
6A	Y002223	16507..(L)	1		Medium R	1627504	1
	Y002224	16508..(R)	1	22	---	1001144	1
6B	Y021872	20183..(L)	1	23	M6x20 DIN933/8.8	1001430	1
	Y021873	20184..(R)	1	24*	M6 type S DIN17222	1912027	1
7	Y002226	16509..	1	25	∅3.9x9.5	1001530	5
8	Y012159	16722..	1	26	∅16	1001598	1
9	Y012670	1672404	1	27	8mm	1001637	2
10	Y012695	1672300	1	28	M6 DIN985	1001680	1
11	X001117	1009204 (L)	1	29*	1/4"-20x1-1/2"	1912025	1
	X001118	1009205 (R)	1	30	X005144	1002085	1
12	X001175	14581..(L)	1	31	Large	1002158	1
	X001176	14582..(R)	1		Medium	1002163	1
13	X001424	1009752	1	32A	M6x35	1006057	3
14A	X001474	1005641	1	32B	M6x35	1006057	2
14B	X005557	2018500	1	33	∅20x∅10.2x1.1	1229600	1
15	X003889	1901094	1	34	M6x20	1006999	2
16	X001470	1397900	1	35	M8x10 DIN916/45H	1008942	1
17	X001476	1008444	2	36A	M6x12mm	1008492	1
18	X003590	1003917	1	36B	M6x12mm	1008492	2
19	X003890	1901095	1	37A	---	1991341 (L)	1
20	M5 DIN934	1001669	1		---	1991342 (R)	1
				37B	X005556	1903947	0.5
				38	M5x40	1902842	1
				39	∅22	1990648	1

Option B complete  
E001017 = 1912114

Option A: Knee protector PU

Option B: Knee protector V40

\* Only for option B

LEG REST BZ7 TYPE 708D				Plan No. / Plannr.
BEENSTEUN BZ7 TYPE 708D				
Scale / Schaal: 3:10	Drafted / Getekend	Verified / Nagezien	Confirmed / Bevestigd	
Date / Datum	13/07/2007	11/09/2017	11/09/2017	
Init. / Init.	SD	SD	PJ	