

* - used only with SE19

EP00488 / C


| | | | | |
|-----------------------|------------|------------|------------|---|
| SE53 ARMREST / TURIOS | | | | Plan No. EP00488 / C |
| Scale: 1:4 | Drafted | Verified | Confirmed |  |
| Date | 12/02/2020 | 30/04/2020 | 30/04/2020 | |
| Initials | DR | DB | TK | |

| Material | | | |
|----------|----------------|-------------|----------|
| No. | Plan | Article | Quantity |
| 1 | - | 2606454 | 1 |
| 2 | Y003626 | 1476554 | 1 |
| 3 | Y009175 | 1641654 | 3 |
| 4 | Y009176 | 1641754 | 1 |
| 5 | YP07839 | 2783154 | 2 |
| 6 | YP08051 | 2605154 | 1 |
| 7A | YP08058 | 2605454 (R) | 1 |
| | YP08059 | 2605554 (L) | 1 |
| 7B* | YP08060 | 2605654 (R) | 1 |
| | YP08061 | 2605754 (L) | 1 |
| 8* | YP08063 | 2605854 | 1 |
| 9 | X005517 | 2014999 | 1 |
| 10 | X005519 | 1911778 | 1 |
| 11 | X005561 | 2021699 | 1 |
| 12 | X005565 | 1912287 | 1 |
| 13 | X005567 | 1912286 (R) | 1 |
| | X005566 | 1912285 (L) | 1 |
| 14 | XP00227 | 2990147 | 2 |
| 15 | Sticker | T001002 | 1 |
| 16 | Washer M6 | 1000973 | 8 |
| 17* | M5x12 | 1001489 | 2 |
| 18 | Tube Cap 20x20 | 1001612 | 2 |
| 19 | Nut M10 | 1001672 | 1 |
| 20 | Nut M5 | 1001679 | 4 |
| 21 | Nut M6 | 1001680 | 6 |
| 22 | M6x40 | 1003111 | 2 |
| 23 | Star Knob | 1003658 | 1 |
| 24 | M5x20 | 1003907 | 4 |
| 25 | M6x40 | 1006111 | 4 |
| 26 | M5x16 | 3991136 | 2 |
| 27 | M5x16 | 1007198 | 2 |

| Material | | | |
|----------|------------------|---------|----------|
| No. | Plan | Article | Quantity |
| 28 | M6x16 | 1008459 | 1 |
| 29 | Nut M6 | 1008492 | 3 |
| 30 | Washer M5 | 1906561 | 2 |
| 31* | M5X5 | 1906562 | 2 |
| 32 | Linx DLX REM 215 | 1910700 | 1 |
| | Linx DLX REM 050 | 1909020 | 1 |
| | Linx SE19 SE39 | 1912296 | 1 |
| | Linx SE39 | 1912301 | 1 |
| 33 | M5x16 | 1907493 | 2 |
| 34* | Light controller | 1915936 | 1 |
| 35 | Spring plunger | 1991432 | 1 |
| 36 | Rivet Ø4 | 1992805 | 2 |
| 37* | M3x8 | 2990309 | 2 |
| 38 | M6x45 | 2994058 | 2 |
| 39 | Tube Cap 16x16 | 2994195 | 2 |

* - used only with SE19

EP00488 / C

| | | | | |
|-----------------------|------------|------------|------------|---|
| SE53 ARMREST / TURIOS | | | | Plan No. |
| | | | | EP00488/ C |
| Scale: 1:4 | Drafted | Verified | Confirmed |  |
| Date | 12/02/2020 | 30/04/2020 | 30/04/2020 | |
| Initials | DR | DB | TK | |