



Material			
No.	Plan	Article	Quantity
1	---	1003299	4 x 0,045 m
2	YP04966	2894904	2
3	XP01322	GP06300	2
4	XP01323	GP06400	2
5	---	---	---
6	D17XD6.4X3 DIN7349	1003321	4
7	M6 DIN127	1004414	4
8	M6X12 ISO7380	1007000	4
9	---	1994568	4
10	---	1003049	4

EP00186 / B

HOME TRAINER TRICYCLES TRENAZER ROWEROWY				Plan No. EP00186/ B
Scale: 1:4	Drafted	Verified	Confirmed	
Date	30/09/2014	12/09/2019	12/09/2019	
Initials	WM	AZ	TK	